



## ITEM TKL-079

<b>Product Name</b>	Digital TENS (Transcutaneous Electrical Nerve Stimulator)	
<b>Package</b>	36PCS/CTN	
<b>Cuft</b>	1.6	
<b>N.W.</b>	9KGS	
<b>G.W.</b>	10KGS	
<b>Description</b>	<ul style="list-style-type: none"> <li>* Relieve soreness in the shoulder</li> <li>* Relieve the paralysis of the peripheral nervous</li> <li>* Massage effect</li> </ul>	
<b>LCD Display</b>	 <div style="display: flex; justify-content: space-around; margin-top: 5px;"> <div style="text-align: left;"> <p>← Output electrode indication</p> <p>← Intensity indication</p> <p>← Time setting: 15/30/60 minutes</p> </div> <div style="text-align: center;"> <p>↑ Mode selection</p> </div> </div>	
<b>Specifications</b>	<ul style="list-style-type: none"> <li>* Dimension: 57mm(W)x118mm(H)x21mm(D)</li> <li>* Output frequency: under 180Hz</li> <li>* Polarity: Automatic change.</li> <li>* Display: Numerical LCD (liquid crystal display)</li> <li>* Button: Rubber made in double colors.</li> <li>* Circuit: Patented PWM (Pulse-Width-Modulation) for DC/DC circuit with high efficiency and resolution.</li> <li>* Protection: Patented ADC (Analogical-Digital-Converter) type of protective circuit for output.</li> <li>* Intensity: 8 power levels selection</li> <li>* Mode: 8 modes selection A/B/C/D/E/F/G/H included multi-programs with intensity in each mode.</li> <li>* Timer: Triple selection 15/30/60 minutes. Auto shut off.</li> <li>* Power supply: Three 1.5V (LR-03/AAA) batteries.</li> <li>* Voltage: 4.5V battery.</li> <li>* Weight: 120 g (battery included)</li> <li>* Ambient Temperature: Operating: 10°C -40°C ; Storage: -20°C -50°C</li> <li>* Humidity: Operating: 20% -95% ; Storage: 10% -95%</li> </ul>	